



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Static pipette holder 340 for 2 Socorex microliter pipettes

€87.98
plus VAT &
Shipping

Product Images

Description

Static pipette holder 340

Stand for micro, macro or multichannel pipettes from Socorex.

Space-saving pipette stand with innovative design for 2 multichannel or single-channel micropipettes.

Holds all Acura models and allows user-friendly access to each instrument.

Heavy duty stand plate with anti-slip feet ensures maximum stability.

Comfortable access to any instrument possible.

Made of robust polyamide, easy to clean.

Additional Information

| | |
|------------------------|------------------|
| No. | SR-0135 |
| Manufacturer (Brand) | Socorex |
| old neoLab article no. | 7-7402 |
| EAN | 4058072366650 |
| Transport temperature | Room temperature |

