



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Static Pipette Holder 340 for 3 Socorex Microliter Pipettes

**€100.78**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Static pipette holder 340

Stand for micro, macro or multichannel pipettes from Socorex.  
Space-saving pipette stand with innovative design for 3 multichannel or single-channel micropipettes.  
Holds all Acura models and allows user-friendly access to each instrument.  
Heavy duty stand plate with anti-slip feet ensures maximum stability.  
Comfortable access to any instrument possible.  
Made of robust polyamide, easy to clean.

## Additional Information

---

No.	SR-0136
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7403
EAN	4058072366667
Transport temperature	Room temperature
Color	gray
Material	Polyamide (PA)
for pipettes / dispensers from	Socorex
Type frame	Pipette rack

