



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Antistatic brush, brush head 20 mm Ø

€23.50
plus VAT &
Shipping

Product Images



Description

Antistatic brushes effortlessly remove dust particles and static charges.
In case of greasiness, simply wash out the fibers and natural hair with soap and lukewarm water.

Additional Information

No.	2-3265
Manufacturer (Brand)	neoLab
EAN	4058072001155
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Form	pin-shaped
Brushes and brushes properties	antistatic
DM Head	20 mm
Strength bristles	soft
Material bristles	Goat hair
Material stem	Polycarbonate (PC)
Type brush	Antistatic brush
Type cleaning tool	Brush. general application
suitable for cleaning of (e.g.)	variable

