



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine cleaning brush for sedimentation vessels, 70 cm long

€35.70
plus VAT &
Shipping

Product Images



Description

Conical cleaning brush made of horsehair/nylon with brush tip for large conical vessels, especially for sedimentation vessels according to Imhoff.

Galvanized steel wire with angled, rotating wooden handle for turning the brush in sedimentation vessels.

The soft bristles ensure fast, thorough cleaning of the inner walls even of heavily soiled sedimentation vessels.

Top diameter: 130 mm

Total length: approx. 70 cm

Brush length: approx. 50 cm

Additional Information

No.	RL-0027
Manufacturer (Brand)	neoLabLine
old neoLab article no.	2-3285
EAN	4058072980658
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	brown
Form	not specified
Length	700 mm
Brushes and brushes properties	flexible
DM Head	130 mm
Strength bristles	soft
Material bristles	Horsehair Polyamide 6.6 (PA 6.6. Nylon)
suitable for cleaning of (e.g.)	Sedimentation Vessels
Length brush	500 mm
Material stem	Wood steel galvanized
Type brush	Brush for special laboratory vessels
Type cleaning tool	Brush. especially for the laboratory

