



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex nylon brush for 0.3, 0.5 and 1 mL glass cylinders for Dosystm self-filling syringes

€34.74
plus VAT &
Shipping

Product Images



Description

Nylon brushes for cylinders

Purpose-designed brushes with head bundle and plastic handle provide efficient cleaning of glass cylinders and tubes.

Various brush sizes available for different cylinder volumes.

Additional Information

No.	SR-0286
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7686
EAN	4058072368043
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Blue White
Form	straight
Material bristles	Polyamide 6.6 (PA 6.6. Nylon)
suitable for cleaning of (e.g.)	Cylinders
Material stem	Steel galvanized
Type brush	Brush for special laboratory vessels
Type cleaning tool	Brush. especially for the laboratory

