



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex nylon brush set for all sizes of glass cylinders, for Dosys™ self-filling syringes

**€34.74**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Nylon brushes for cylinders

Purpose-designed brushes with head bundle and plastic handle provide efficient cleaning of glass cylinders and tubes.

Various brush sizes available for different cylinder volumes.

## Additional Information

---

No.	SR-0289
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7689
EAN	4058072368067
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Blue White
Form	straight
Brushes and brushes properties	in set
Material bristles	Polyamide 6.6 (PA 6.6. Nylon)
suitable for cleaning of (e.g.)	Cylinders
Material stem	Steel galvanized
Type brush	Brush for special laboratory vessels
Type cleaning tool	Brush. especially for the laboratory

