




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Articulated Glass Brush

**€14.90**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Round piston cleaning brush, universally usable.

The joint on the brush also allows the upper inner neck of the piston to be reached.

The tiltable brush body allows every part of the inner wall of the piston to be reached and cleaned.

Natural bristles adapted to the curve of the inner wall of the piston ensure thorough cleaning even with heavy soiling.

The brush fits in every piston with ground joint NS 29/32 or larger.

The bristles used are boiled at 100°C during production.

Good grip due to the sturdy, long wooden handle.

The joint glass brush is also well suited for pre-cleaning.

## Additional Information

---

No.	1-6010
Manufacturer (Brand)	neoLabLine
EAN	4058072036140
autoclavable at 121°C	No
Color	wood color
Wide	22 mm
Form	not specified
Height	20 mm
Brushes and brushes properties	in set
Strength bristles	hard
Material bristles	Pig bristles
suitable for cleaning of (e.g.)	Test tubes flasks cylinders variable
Spare part	1-6011
Length bristles	10 mm
Length brush	45 mm
Material stem	Wood
Type brush	Brush for special laboratory vessels
Type cleaning tool	Brush. general application

