



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY eye emergency station, wall box 265 x 227 x 110 mm

€75.12
plus VAT &
Shipping

Product Images



Description

Eye emergency stations for use in workplaces where chemicals and also foreign bodies can get into the eye. For wall mounting, including 200 ml pH Neutral rinsing solution for chemical substances and 500 ml eye rinsing solution for foreign bodies.

Even after rinsing with pH Neutral, it is recommended to continue rinsing with the eye rinsing solution until medical treatment is received.

Each with integrated pictogram and mirror, including eye emergency plan.

Rinsing time approx. 7 minutes each.

Dust-tight, closed wall box

Additional Information

No.	FG-0043
Manufacturer (Brand)	B-SAFETY
old neoLab article no.	1-3034
EAN	4058072254162
Transport temperature	Room temperature
Color	green
Volume	700 ml
sterile	No
Wide	110 mm
Height	265 mm
Length	227 mm
Depth	227 mm
Showers and eye wash properties	suitable for wall mounting
Durability	3 years
Type solution	Ophthalmic rinsing solutions (foreign bodies) Ophthalmic rinsing solutions (chemical substances)
filled according to	DIN 12930
Contents eye rinsing solution	Rinsing solution for foreign bodies pH-Neutral Rinsing solution
Wall suspension	closed
standardized according to DIN standard (DIN EN)	DIN EN 15154-4:2010-04

