



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® First Aid Kit, filled, according to DIN 13157

€155.20
plus VAT &
Shipping

Product Images



Description

In addition to the standard contents according to DIN 13157, this first aid kit contains valuable supplements such as ^{aluderm}® special dressings, special pressure dressing packs, a hand cleaning set and extended plaster equipment, and much more.

The sterile dressing materials are CE-marked, with the date of manufacture and a shelf life of 20 years.

Convenient wall mount with 90° stop lock. Can also be opened in the wall holder. Handy grip keys. Transparent cover plates, circumferential rubber seal, adjustable interior compartments, smooth-running twist locks, reinforced carrying handle, sealing device, security seal.

ABS plastic, regenerate-free, cadmium-free coloring, approved for medical devices, pharmaceuticals and food. Dimensionally stable, impact-resistant, temperature-resistant, splash-proof, dust-proof, maintenance-free.

Dimensions: 310 x 210 x 130 mm (L x W x D)

Additional Information

No.	2-2100
Manufacturer (Brand)	neoLab
EAN	4058072059774
Transport temperature	Room temperature
Color	orange
Material	Acrylonitrile butadiene styrene (ABS)
Wide	130 mm
Height	210 mm
Length	310 mm
Depth	310 mm
Type equipment	Standard
Type first aid equipment	First aid kit
filled according to	DIN 13157

