



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## B-SAFETY PPE cabinet CLASSIC - for storage of PPE

**€172.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

B-SAFETY PPE cabinet CLASSIC is the ideal solution to store personal protective equipment. The cabinet is not only practical, spacious and reliable, but also designed for safety in the workplace. B-SAFETY PPE cabinet CLASSIC is made of first-class materials and excellent quality. It is injection molded from polypropylene and guarantees excellent durability and high reliability in any situation. PPE cabinet CLASSIC is the right choice for those who are looking for a practical and safe product. Suitable for equipping with various PPE items from hard hats to safety goggles. Product features - good resistance to hydrocarbons and lubricating oils, good UV resistance - includes two shelves for sorted storage of various items - two inspection windows with perimeter seal for good water resistance - internal seal to protect against dust and water penetration - includes two ergonomic nylon latches for good closure - Material: Polypropylene, nylon, polycarbonate - Color: green/black - Mounting type: wall-mounted - Protection type: IP59K - Weight: 7.2 kg - Dimensions (H x W x D): 735 x 616 x 270 mm

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | FG-0250          |
| Manufacturer (Brand)  | B-SAFETY         |
| EAN                   | 4066292559028    |
| Transport temperature | Room temperature |

