



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® dressing cabinet with filling, 290 x 360 x 140 mm

€130.50
plus VAT &
Shipping

Product Images



Description

Practical first-aid cabinet, very solidly manufactured.
Due to its compact design, it is particularly suitable for smaller operating units.
Made of white painted sheet steel with lockable door.

Filling DIN 13157 according to ASR A4.3
Dimensions: 302 x 362 x 140 mm (W x H x D)

Additional Information

No.	2-2102
Manufacturer (Brand)	neoLab
EAN	4058072059798
Transport temperature	Room temperature
Color	White
Material	Steel powder coated
Wide	140 mm
Height	302 mm
Length	362 mm
Depth	362 mm
filled according to	DIN 13157 according to ASR A4.3
Type equipment	Standard
Type first aid equipment	Dressing cabinet

