



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY plaster dispenser system with 1x45 waterproof and 1x45 textile plasters

€28.68
plus VAT &
Shipping

Product Images



Description

The QuickFix wound plaster dispenser provides quick and hygienic access to plasters directly at the workplace - simply pull and the plaster is immediately ready for use. Just apply the plaster to the wound and it's back to work. The QuickFix plaster dispenser system is supplied with plasters, fastening material and key for theft protection

2-0311: Plaster dispenser system with 1x45 waterproof and 1x45 textile plasters

Additional Information

No.	FG-0072
Manufacturer (Brand)	B-SAFETY
old neoLab article no.	2-0311
EAN	4058072254193
Transport temperature	Room temperature
Color	beige
Number of parts in set	45
Refill package/article	FG-0071 FG-0074

