



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic Resin Spray, 200 ml

€17.10
plus VAT &
Shipping



Description

Protective film for printed circuits, safely protects against leakage currents and fine connections.

Reliable protection against oxidation and moisture protected.

Insulates cables and wires.

The highly insulating, crystal clear acrylic resin lacquer is resistant to acids, alkalis, alcohol, atmospheric influences and extreme temperatures.

Additional Information

No.	1-6238
Manufacturer (Brand)	neoLab
EAN	4058072043124
Volume	200 ml
Type special spray	Plastic spray

