



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Labosol S liquid, bottle 1 l

€32.50
plus VAT &
Shipping

Product Images



Description

The ingredients of Labosol® S are constantly monitored.

Intensive cleaning concentrate for soiled laboratory glassware, porcelain, instruments and material made of plastic and rubber.

Cleans thoroughly and automatically without attacking the material to be cleaned or the skin. Can also be used for ultrasound.

0.5-2% solution depending on contamination. Cleaning power can be increased by heating the bath.

Substitute for chromosulfuric acid, without chlorine and caustic alkalis, biodegradable.

Additional Information

No.	1-6091
Manufacturer (Brand)	neoLabLine
EAN	4058072036577
Volume	1000 ml
Detergent properties	biodegradable
Medium/ Consistency	Liquid
Special cleaner for cleaning (objects/ devices)	Ultrasonic equipment Laboratory equipment Laboratory instruments
Type detergent	Surface cleaner
bottled in	Bottle with screw cap

