



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Labosol U, ultrasonic cleaner, canister 5 l

€175.80
plus VAT &
Shipping

Product Images



Description

The Labosol® ultrasonic^{bath} cleaning concentrate intensifies the cavitation of the ultrasound. Therefore cleans all surfaces quickly, powerfully and gently from oil, grease, wax, resin, silicone, blood, protein and others.

Labosol® U is water soluble, non-flammable and has a fast degassing effect.

Very low application concentrations, can be used several times. Only weakly foaming.

Can replace the dangerous and difficult to dispose of chromosulfuric acid in most cases.

Additional Information

No.	1-6788
Manufacturer (Brand)	neoLabLine
EAN	4058072044367
Volume	5000 ml
Medium/ Consistency	Liquid
Type detergent	Special cleaner for special objects and equipment
bottled in	Bottle with screw cap

