



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Labothermol® K -20 to +150°C, canister 5 l

€156.60 As low  
as  
€137.60  
plus VAT &  
Shipping

### Product Images



## Description

---

Safety refrigeration and heat transfer fluid.

Water soluble bath fluid for refrigeration and heat thermostats.

- washable with water
- facilitates cleaning of equipment
- water that accidentally enters the heating bath does not splash, even at temperatures above 100°C, because it can be safely separated.

Modified, polyvalent, aliphatic alcohol mixed with 15% water, almost colorless.

Flash point according to DIN 51 758 above 255°C.

Temperature range: -20 °C to +150 °C

Pumpable down **to 0 °C**

## Additional Information

---

No.	1-6065
Manufacturer (Brand)	neoLabLine
EAN	4058072036461
Temp MAX	150 °C
Temp MIN	-20 °C
Volume	5000 ml

