



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**neolab® ESD cleaning sticks, pointed-conical,  
head PU 10x3.3 mm Ø, handle PP 65 mm, total  
length 75 mm**

**€53.80  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

The **antistatic ESD cleaning sticks**, or **ESD Swaps**, are suitable for cleaning surfaces and poorly accessible areas. Their surface resistance is between  $1 \times 10^6 \Omega \leq R_s \leq 1 \times 10^{10} \Omega$ . The heads of the rods are made of polyurethane foam PU (100 ppi) or polyester. The cleaning head is attached to the stem without adhesive. For product cleaning of lenses, sensors, semiconductors and microchips, can also be used in clean rooms.

- ideal for use with numerous cleaning agents and solvents

## Additional Information

---

No.	2-1059
Manufacturer (Brand)	neoLabLine
EAN	4058072055325
Color	White
Material	Polyurethane (PU)
Length	75 mm
Figure	3
DM Head	3.3 mm
Shape lace	pointed-conical
Type cleaning tool	Micro / pipe cleaner

