



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® silicone gasket with PTFE cuff bore 26 mm No.2-3053

€182.90
plus VAT &
Shipping

Product Images



Description

Silicone seal for screw cap GL 45 with bore, red

With vulcanized PTFE cuff.

Bore: 26 mm

Accessories for screw cap (Order No. 2-3051)

Additional Information

No.	SD-2745
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3052
EAN	4058072363529
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	white transparent (yellow)
Material	Silicone Polytetrafluoroethylene (PTFE)
DM Hole	26 mm
DM top	42 mm
DM below	42 mm
Material gasket / septum	Silicone
Material cuff	Polytetrafluoroethylene (PTFE)

