



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® silicone gaskets (VMQ) for piercing (Septa) for screw joint closure, GL 14

€90.00
plus VAT &
Shipping

Product Images



Description

DURAN® GL 14 silicone septum seal, peroxide cured silicone, VMQ Fits screw joint closure made of PBT. Thermal resistance: 130 °C (steam) and 200 °C (heat). The silicone is peroxide cured. Example applications: Injection or withdrawal of media. Special features and benefits: Fits PBT bore screw cap Compatible with GL bottle threads Peroxide cross-linked silicone Application: Aseptic reagent delivery or culture sampling using syringe and needle, Aseptic reagent delivery or culture sampling using syringe and needle

Additional Information

| | |
|-----------------------|------------------|
| No. | SD-2775 |
| Manufacturer (Brand) | DURAN |
| VGKL number | 132924605 |
| EAN | 4058072363772 |
| Manufacturer EAN | 4032051064058 |
| Gross weight | 0.001 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| Net weight | 0.00100 kg |
| Material | Silicone |
| Country of origin | Germany |
| sterile | No |
| Region of origin | Lower Saxony |
| Customs tariff number | 39269099 |
| DM below | 12 mm |

