



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Septum ND24 silicone natural/PTFE beige, 1000 pcs/pack

€503.90
plus VAT &
Shipping

Product Images



Description

Unmounted septa ND24.

Septum Silicone/PTFE

Temperature resistant from -60°C to 200°C

For easy penetration

For critical analysis

Thickness 3.2 mm

Shore hardness 45° shore A

Additional Information

No.	7-0866
Manufacturer (Brand)	neochrom
EAN	4058072037345
Transport temperature	Room temperature
Thickness septum	3.2 mm
Hardness	45 shore A
TBST MAX closure/ lid	200 °C
TBST MIN Septum	-60 °C

