



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Septum Silicone/PTFE, 8 mm Ø, 100 pcs./pack

€47.00
plus VAT &
Shipping

Product Images



Description

Unassembled septa for ND8 threaded bottles.

Additional Information

No.	2-7204
Manufacturer (Brand)	neochrom
EAN	4058072009014
Transport temperature	Room temperature
Thickness septum	1.3 mm
Hardness	45 shore A
TBST MAX closure/ lid	200 °C
TBST MIN Septum	-60 °C

