



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Septum ND8 silicone white/PTFE blue, slotted, 1000 pcs./pack

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Unassembled septa for ND8 threaded bottles.

Additional Information

No.	7-0670
Manufacturer (Brand)	neochrom
EAN	4058072029340
Transport temperature	Room temperature
Hardness	55 shore A
Thickness septum	0.9 mm
TBST MAX closure/ lid	200 °C
TBST MIN Septum	-60 °C

