



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neochrom® Septum ND20, silicone blue  
transparent /PTFE transparent, 1000 pcs./pack**

**€513.00  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

Septa for use with ND20 flare caps.

Septum Silicone/PTFE

Temperature resistant from -60°C to 200°C

For easy penetration

For critical analyses

## Additional Information

---

No.	7-0804
Manufacturer (Brand)	neochrom
EAN	4058072030506
Transport temperature	Room temperature
Hardness	45 shore A
Thickness septum	3 mm
TBST MAX closure/ lid	200 °C
TBST MIN Septum	-60 °C

