



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz three-way plug valve with PTFE plug, valve NS 14.5, bore 1.5 mm

**€45.25**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Three-way Plug Valves with PTFE Plug  
According to DIN 12554 and in addition to DIN.  
T-bore with PTFE plug and threaded valve lock.  
Tightness according to DIN 12540-1.

## Additional Information

---

No.	LZ-0377
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2376
EAN	4058072301194
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	1.5 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44695
Type cock / valve	Tap 3-way
Type connection	Glass spigot

