



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz three-way plug valve with PTFE plug, valve NS 14.5, bore 2.5 mm

€46.88
plus VAT &
Shipping

Product Images



Description

Three-way Plug Valves with PTFE Plug
According to DIN 12554 and in addition to DIN.
T-bore with PTFE plug and threaded valve lock.
Tightness according to DIN 12540-1.

Additional Information

No.	LZ-0378
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2377
EAN	4058072301200
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44695
Type cock / valve	Tap 3-way
Type connection	Glass spigot

