



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz three-way plug valve with solid glass plug, valve NS 18.8, bore 2.5 mm

€44.19
plus VAT &
Shipping

Product Images



Description

Three-way plug valves with glass plugs

According to DIN 12554 and in addition to DIN.

T-bore with solid glass plug and threaded valve lock.

Tightness according to DIN 12540-1.

Additional Information

No.	LZ-0373
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2372
EAN	4058072301156
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	44791
Type cock / valve	Tap 3-way
Type connection	Glass spigot

