



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz three-way plug valve acc. to Czako with PTFE plug, valve NS 14.5, bore 2.5 mm

**€85.63**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Three-way Plug Valves acc. to Czako with PTFE Plug

According to DIN 12563; angular bore: 120°.

With threaded valve lock; tightness according to DIN 12540-1.

## Additional Information

---

No.	LZ-0384
Manufacturer (Brand)	Lenz
old neolab article no.	B-2384
EAN	4058072301262
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44695
Type cock / valve	Tap 3-way
Type connection	Glass spigot

