



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene two-way stopcock (PP), with 4 mm hole, TEF chicks

€51.00 As low as
€41.33
plus VAT &
Shipping

Product Images



Description

Two-way faucet

The two ends of the tube are shaped as hose olives.

Ideal as a replacement for glass faucets, unbreakable, extremely tight, do not need to be greased.

Additional Information

No.	CD-1777
Manufacturer (Brand)	Nalgene
old neolab article no.	1-0046
EAN	4058072219529
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky) white black
sterile	No
Roosters properties	with special chemical resistance
Number of hose olives	3
Bore	4 mm
Material tap	Polypropylene (PP)
Material chicks	Polytetrafluoroethylene (PTFE. Teflon)
Type cock / valve	Tap 2-way
for DM hose inside 1	6.0-9.0 mm
Type connection	Hose olive

