



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene three-way stopcock (PP), with 4 mm bore, TEF plug

**€59.70** As low  
as  
**€48.25**  
plus VAT &  
Shipping

### Product Images



## Description

---

Three-way faucet

The three pipe ends are shaped as hose olives.

Ideal as a replacement for glass faucets, unbreakable, extremely tight, do not need to be greased.

## Additional Information

---

|                        |  |
|------------------------|--|
| No.                    | CD-1779                                |
| Manufacturer (Brand)   | Nalgene                                |
| old neoLab article no. | 1-0044                                 |
| EAN                    | 4058072219536                          |
| Transport temperature  | Room temperature                       |
| autoclavable at 121°C  | Yes                                    |
| Color                  | transparent (milky) white black        |
| sterile                | No                                     |
| Roosters properties    | with special chemical resistance       |
| Number of hose olives  | 2                                      |
| Bore                   | 4 mm                                   |
| Material tap           | Polypropylene (PP)                     |
| Material chicks        | Polytetrafluoroethylene (PTFE. Teflon) |
| Type cock / valve      | Tap 3-way                              |
| for DM hose inside 1   | 6.0-9.0 mm                             |
| Type connection        | Hose olive                             |

