



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® metal three-way stopcock with LL connections, 2 LLW, 1 LLM

€55.40
plus VAT &
Shipping

Product Images



Description

Metal stopcock with Luer-Lock connections.

The stopcock plug is held in place in the stopcock bed by means of a spring.

Silicone grease is used as a lubricant.

Additional Information

No.	2-6286
Manufacturer (Brand)	neoLab
EAN	4058072007942
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
sterile	No
Roosters properties	with special chemical resistance with special pressure resistance
Quantity Luer lock male	1
Quantity Luer lock female	2
Bore	1.8 mm
Material tap	Brass chrome plated
Material chicks	Brass chrome plated
Type cock / valve	Tap 3-way
Type connection	Luer lock female Luer lock male

