



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Mini three-way stopcock, 2 x Luer female, 1 x Luer male

€3.81
plus VAT &
Shipping

Product Images



Description

Small three-way stopcock made of plastic, sterile, for connecting syringes or cannulas.

Additional Information

No.	2-4190
Manufacturer (Brand)	neoLab
EAN	4058072002404
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
sterile	Yes
Quantity Luer cone female	2
Quantity Luer lock male	1
Material tap	High-density polyethylene (HDPE)
Material chicks	High-density polyethylene (HDPE)
Type cock / valve	Tap 3-way
Type connection	Luer cone female Luer lock male

