



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz Straight burette stopcock NS 12.5, bore 1.5 mm, with solid glass plug

€23.85
plus VAT &
Shipping

Product Images



Description

Burette taps with faceted spout tip.

Additional Information

No.	LZ-0386
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2386
EAN	4058072301286
Transport temperature	Room temperature
Material tap	Borosilicate glass 3.3
Material chicks	Borosilicate glass 3.3
Standard ground	44693
Type Connection Output	Glass spigot
Type Connection Input	Glass spigot

