



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz Lateral burette stopcock NS 12.5, bore 1.5 mm, with PTFE plug

**€33.18**  
plus VAT &  
Shipping

### Product Images



## Description

---

**Burette taps** with faceted spout tip.

## Additional Information

---

No.	LZ-0390
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2390
EAN	4058072301323
Transport temperature	Room temperature
Material tap	Borosilicate glass 3.3
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	44693
Type Connection Output	Glass spigot
Type Connection Input	Glass spigot

