



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz vacuum plug valve, one-way, straight, NS 14.5

€48.37
plus VAT &
Shipping

Product Images



Description

One-way **Vacuum Plug Valves**, Straight (Bag Valves)

Straight-through stopcocks with straight neck tubes; with hollow glass plug.

Additional Information

| | |
|------------------------|---------------------|
| No. | LZ-0396 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-2396 |
| EAN | 4058072301385 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| sterile | No |
| Bore | 2.5 mm |
| Material tap | Borosilicate glass |
| Material chicks | Borosilicate glass |
| Standard ground | 44695 |
| Type cock / valve | Tap 1-Way |
| Type connection | Glass spigot |

