



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz vacuum plug valve, one-way, right-angled, valve NS 21.5

**€60.43**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

One-way **vacuum plug valves** (angle valves) with right-angled attachment tubes and hollow glass plug.

## Additional Information

---

No.	LZ-0403
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2403
EAN	4058072301453
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	6 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	44702
Type cock / valve	Tap 1-Way
Type connection	Glass spigot

