



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz one-way plug valve with PTFE plug, valve NS 12.5, bore 2.5 mm

**€283.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### One-way Plug Valves with PTFE Plug

According to DIN 12541-1 and in addition to the DIN standard. With PTFE plug and threaded valve lock.

One straight bore; material: borosilicate glass 3.3.

The tightness meets the requirements of DIN 12540-1.

## Additional Information

---

No.	LZ-0339
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2338
EAN	4058072300814
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	44693
Type connection	Glass spigot
Type cock / valve	Tap 1-Way

