



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz one-way plug valve with PTFE plug, valve NS 18.8, bore 4.0 mm

€43.59
plus VAT &
Shipping

Product Images



Description

One-way Plug Valves with PTFE Plug

According to DIN 12541-1 and in addition to the DIN standard. With PTFE plug and threaded valve lock.

One straight bore; material: borosilicate glass 3.3.

The tightness meets the requirements of DIN 12540-1.

Additional Information

No.	LZ-0342
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2341
EAN	4058072300845
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	4 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44791
Type cock / valve	Tap 1-Way
Type connection	Glass spigot

