



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PP replacement stopcock for desiccator No. 2-5311/14

€52.50
plus VAT &
Shipping

Product Images



Description

3-way stopcock

Housing made of PP, closure and sealing ring made of PTFE

Spare part for desiccator with stopcock (Order No. 2-5311 and 2-5313)

Additional Information

No.	2-5315
Manufacturer (Brand)	neolab
EAN	4058072004613
Transport temperature	Room temperature
Material tap	Polytetrafluoroethylene (PTFE, Teflon)
Type Connection Output	Hose olive
Type Connection Input	Tube
for DM hose inside l	8 mm
for type tube	Cover tube side tube

