



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex fluid control system with stopcock for Calibrex™ 525/530 dispensers up to 10 ml

€83.72
plus VAT &
Shipping

Product Images



Description

Fluid control system with stopcock

Compatible with any Calibrex™ 525/530.

In addition to the dispense position, this element allows initial filling and recirculation of reagents without their loss. The locking position of the stopcock provides additional safety.

Completely autoclavable.

Additional Information

No.	SR-0237
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7614
EAN	4058072367558
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Grey Black
sterile	Yes
Bore	10 mm
Material tap	Polypropylene (PP)
Material chucks	Polytetrafluoroethylene (PTFE)
Type cock / valve	Tap 1-Way
Type connection	Screw adapter

