



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Fluid Control System with Stopcock for Calibrex™ 525/530 Dispensers 25,50,100ml

**€83.72**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Fluid control system with stopcock

Compatible with any Calibrex™ 525/530.

In addition to the dispense position, this element allows initial filling and recirculation of reagents without their loss. The locking position of the stopcock provides additional safety.

Completely autoclavable.

## Additional Information

---

No.	SR-0238
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7615
EAN	4058072367565
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Grey Black
sterile	Yes
Bore	15 mm
Material tap	Polypropylene (PP)
Material chucks	Polytetrafluoroethylene (PTFE)
Type cock / valve	Tap 1-Way
Type connection	Screw adapter

