



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz glass stopcock, solid, NS 12.5, bore 2.5 mm, SBW

€9.75
plus VAT &
Shipping



Description

Solid **glass stopcock** with thread lock.
For one-way plug valves; made of glass.

Additional Information

No.	LZ-0330
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2329
EAN	4058072300722
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Bore	2.5 mm
Material chucks	Borosilicate glass
Standard ground	44693

