




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz PTFE stopcock NS 12.5, bore 1.5 mm

**€11.20**  
**plus VAT &**  
**Shipping**



## Description

---

PTFE stopcock plug with thread and threaded stopcock retainer.  
For one-way plug valves.

## Additional Information

---

No.	LZ-0347
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2346
EAN	4058072300890
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
sterile	No
Bore	1.5 mm
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	44693

