



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz PTFE stopcock NS 12.5, bore 2.5 mm

€11.20
plus VAT &
Shipping



Description

PTFE stopcock plug with thread and threaded stopcock retainer.
For one-way plug valves.

Additional Information

No.	LZ-0348
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2347
EAN	4058072300906
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
sterile	No
Bore	2.5 mm
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	44693

