



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Braun Petzold flow-through tap 1x Luer lock female/male each

€1.20
plus VAT &
Shipping

Product Images



Description

Plastic valves
with end caps, rotatable LL union nut, pressure tested.

Additional Information

No.	TZ-2762
Manufacturer (Brand)	Braun Petzold
old neoLab article no.	2-6358
EAN	4058072398101
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
sterile	Yes
Figure	1
Quantity Luer lock male	1
Quantity Luer lock female	1
Material tap	Polycarbonate (PC)
Material chicks	High-density polyethylene (HDPE)
Type connection	Luer lock female Luer lock male
Type cock / valve	Tap 2-way

