



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz Straight valve stopcock made of PTFE, valve opening 0 - 2.5 mm

Special Price

€15.10 was

€32.47

**plus VAT &
Shipping**



Description

Valve cocks, straight

Smooth-running valve stem made of PTFE,

Valve opening: 0-2,5 mm, with safety device against unintentional unscrewing.

Screw cap and safety device made of PP.

Service temperature: -20 to +150 °C

Additional Information

No.	LZ-0360
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2359
EAN	4058072301026
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass DURAN
Material chucks	Polytetrafluoroethylene (PTFE)
Type connection	Glass spigot
Type cock / valve	Valve tap

