




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz replacement stopcock NS 24 for desiccator with lid tube

€50.09
plus VAT &
Shipping

Product Images



Description

Desiccator replacement valves

With core NS 24, Duran, PTFE spindle

Order No. E-1105: for cover tube

Additional Information

No.	LZ-0394
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2394
EAN	4058072301361
Transport temperature	Room temperature
Material tap	Borosilicate glass
Material valve stem	Polytetrafluoroethylene (PTFE)
Standard ground	24/29
Type Connection Output	Hose olive
Type Connection Input	Standard ground
for DM hose inside 1	8 mm
for type tube	Cover tube

