



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tap bank 2-fold, polysulfone, 3 x LL female, 1 x LL male

Special Price

€4.50 was

€6.75

**plus VAT &
Shipping**

Product Images



Description

Flat design with cross handle

With end caps, pressure-tested, color-coded tap connections.

Color: transparent

Additional Information

No.	2-6355
Manufacturer (Brand)	neoLab
EAN	4058072073756
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
sterile	No
Quantity Luer lock male	1
Quantity Luer lock female	3
Material tap	Polysulfone (PSU)
Material chicks	High-density polyethylene (HDPE)
Type cock / valve	Tap bench 2-fold
Type connection	Luer lock female Luer lock male

