



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz two-way plug valve with solid glass plug, valve NS 14.5, bore 2.5 mm

€59.07
plus VAT &
Shipping

Product Images



Description

Two-way plug valves with solid glass plug

With solid glass plug and threaded cock locking device; two oblique bores.
Tightness meets the requirements of DIN 12540-1.

Additional Information

No.	LZ-0363
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2362
EAN	4058072301057
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	44695
Type cock / valve	Tap 2-way
Type connection	Glass spigot

