



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz two-way plug valve with solid glass plug, valve NS 18.8, bore 4.0 mm

**€67.98**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Two-way plug valves** with solid glass plug

With solid glass plug and threaded cock locking device; two oblique bores.  
Tightness meets the requirements of DIN 12540-1.

## Additional Information

---

No.	LZ-0364
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2363
EAN	4058072301064
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	4 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	44791
Type connection	Glass spigot
Type cock / valve	Tap 2-way

