



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz two-way plug valve with solid glass plug, valve NS 14.5, bore 1.5 mm, capillary attachment tubes

**€72.84**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Two-way plug valves with solid glass plug  
With threaded cock locking device; two oblique bores.  
Tightness meets the requirements of DIN 12540-1.  
With capillary attachment tubes.

## Additional Information

---

No.	LZ-0367
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2366
EAN	4058072301095
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	1.5 mm
Material tap	Borosilicate glass
Material chucks	Borosilicate glass
Standard ground	44695
Type cock / valve	Tap 2-way
Type connection	Glass spigot

